



HIV Antibody Rapid Test

SALIVA samples

BENEFITS

- Rapid test results within 30 minutes
- ELISA-like method
- Complete kit: all materials are included in the box
- Easy, visual interpretation of the results
- Detection of both HIV-1 and HIV-2 antibodies
- Clinical sensitivity: **99.09%**
- Clinical specificity: **99.79%**




New rapid method for detection of HIV 1+2 antibody in oral mucosal transudate fluid


HIV infection can be detected by testing for HIV antibody in serum, urine, oral mucosal transudate fluid, semen, tear of individuals infected by HIV.

This highly sensitive test can effectively and rapidly detect antibody to HIV-1 and HIV-2 in oral mucosal transudate fluid.

Test procedures:




Sample preparation




Adding sample



Adding reagents



Negative



Positive